



Shanghai Joy Crown Industry Co., Ltd.
Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,
Shanghai 200125, P.R. China
Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572
Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Elastic Cohesive Bandage Specifications

Catalogue No. /Reference No. /Products Code	JC04017
Material	60% Non woven, 35% latex or latex free glue, 5% spandex
Weight per Sqm	Fabric Weight / m2: 30gsm; Gluing Weight/m2: 40gsm;
Color	Colored, Printing, Cartoon, Camouflage
Size	Width: 2.5CM, 5CM, 7.5CM, 10CM, 15CM, 20CM, 25CM, 30CM, etc; Length (stretched): 4m, 4.5m, 20m, etc;
Elasticity	210%~220%
Clips	Without Clip or Fastener
Package	Non Sterile: 1 Roll/Candy Bag or Cellophane
Sterile or Not	No Sterile and Sterile
Rules & Classification	Rule 1, Class I



Shanghai Joy Crown Industry Co., Ltd.
Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,
Shanghai 200125, P.R. China
Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572
Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Elastic Cohesive Bandage Specifications


Reference No/ Product Code	Description	Size	Latex or Latex Free	Unit Package	Outside/Cart on Package
JC040171	Elastic Cohesive Bandage, Latex, Tan, Blue, White, Red, Green, Yellow, Black, Pink, etc or OEM printings, 1 Roll/Candy Bag, with/out OEM printing	2.5cm x 4.5m	Latex	48 Rolls/Box	576 Rolls
JC040172		5cm x 4.5m	Latex	24 Rolls/Box	288 Rolls
JC040173		7.5cm x 4.5m	Latex	16 Rolls/Box	192 Rolls
JC040174		10cm x 4.5m	Latex	12 Rolls/Box	144 Rolls
JC040175		15cm x 4.5m	Latex	12 Rolls/Box	96 Rolls
JC040176	Non Woven Cohesive Bandage, Latex Free, Tan, Blue, White, Red, Green, Yellow, Black, Pink, etc or OEM printings, 1 Roll/Candy Bag, with/out OEM printing	2.5cm x 4.5m	Latex Free	48 Rolls/Box	576 Rolls
JC040177		5cm x 4.5m	Latex Free	24 Rolls/Box	288 Rolls
JC040178		7.5cm x 4.5m	Latex Free	16 Rolls/Box	192 Rolls
JC040179		10cm x 4.5m	Latex Free	12 Rolls/Box	144 Rolls
JC0401710		15cm x 4.5m	Latex Free	12 Rolls/Box	96 Rolls

Note: *We also provide following sizes:
Width of 2cm,4cm,6cm,8cm,10cm,12cm, 3.8cm and etc; Length of 2M, 2.7M, 6.9M, 20M and etc. If you have any requirement of other sizes, please contact for more details.

Cohesive Bandage-Pure Color



Snow White-SW01




Beige-BG02



Rose Red-RD03



Olive Green-OG04



Hyacinth Blue-HB05



Black-BK06



Sky Blue-SB07




Lavender Purple-LP08




Neon Green-NG09



Neon Pink-NP10



Navy Blue-NB11



Neon Orange-NO12



Lemon Yellow-LY13




Forest Green-FG14



Coffee Brown-CB15



Peacock Blue-PB16



Jade Green-JG17



Aquamarine-AQ18



Light Purple-LP19



Dark Purple-DP20



Purplish Red-PR21




Orange-OG22






Lake Blue-LK23



Light Brown-LB24



Light Brown-DB25



Cohesive Bandage-Print-Pattern



Captain America-CA31



Super Man-SP32



Batman-BM33




Yellow Smiley-YM34



Red Heart-RH35



Blue Check-BC36



White Skeleton-WS37



Yellow Skeleton-YS38



Black Skeleton-BS39



Green Oval-GO40



Small Green Paw-GP41



Black Horse-GH42




Paw Print-FPD43



Paw Print-WP44



Paw Print-FPD45



Paw Print-RP46



Paw Print-YP47



Paw Print-PPD48



Paw Print-RP49



Paw Print-BPD50




Paw Print-BP51



Paw Print-YPD52



Paw Print-OPD53



Zebra-ZB54

Cohesive Bandage-Camouflage



Debon Camo-DB57



Forest Camo-FR58



Highland Camo-HL59




Sunrise Camo-SU60



Sky Camo-SK61



Ocean Camo-OC62



War Camo-WR63



Jungle Camo-WL64



ACU Camo-AC65




Desert Camo-DS66



Lightgreen Camo-LG67




Snow Camo-SN68



Wild Camo-WL69



Marsh Camo-MR70



Maple Camo-MP71



Grass Camo-GS72

